# Connected Scale How can I set a weight goal?

You can set a weight goal in the ARCHOS Connected Self application. First you will need to press the section dedicated to the Connected Scale.

Page 1 / 5 (c) 2024 ARCHOS <noreply@archos.com> | 2024-04-26 15:30 URL: http://faq.archos.com//index.php?action=artikel&cat=132&id=207&artlang=en

â	ARCHOS Con	nected Self 👤	*
8	David		
Weight		06/06/2014 09:30 76.9 kg	>
	Share	Edit	

In this new screen touch the "goal" button.

Page 2 / 5 (c) 2024 ARCHOS <noreply@archos.com> | 2024-04-26 15:30 URL: http://faq.archos.com//index.php?action=artikel&cat=132&id=207&artlang=en



Page 3 / 5 (c) 2024 ARCHOS <noreply@archos.com> | 2024-04-26 15:30 URL: http://faq.archos.com//index.php?action=artikel&cat=132&id=207&artlang=en

You will need to activate the Goal mode in the top right.



Page 4 / 5 (c) 2024 ARCHOS <noreply@archos.com> | 2024-04-26 15:30 URL: http://faq.archos.com//index.php?action=artikel&cat=132&id=207&artlang=en

You then need to set your target weight goal by using the + and - buttons or drag the ruler to your desired weight then set the date you want to attain your target. Unique solution ID: #1399 Author: Admin FAQ Last update: 2014-07-23 18:05