

# Blood Pressure Monitor

## What do I need on my smart device to use the Blood Pressure Monitor?

In order to be able to sync data from your Blood Pressure Monitor to your Smart Device you will need:

- Either an Android device running Android 4.3 (or later) with Bluetooth® Smart technology, or an iPad® 3 / iPhone® 4S (or later) running iOS 5 to 8.4.1.



- The **ARCHOS Connected Self** application installed - you can download it from [Google Play™](#) and [iTunes®](#)

Unique solution ID: #1401

Author: Admin FAQ

Last update: 2016-01-25 16:07